

## 231/3 BIOLOGY PRACTICAL

### CONFIDENTIAL INSTRUCTIONS TO SCHOOLS

---

*Provide each student with the following;*

- Four (4) dry clean test tubes.
- Test tube holder
- Means of heating.
- Iodine solution.
- Benedicts solution.
- 10% sodium hydroxide solution.
- 1% copper (II) sulphate solution.
- 0.1 % DCPIP solution.
- Measuring cylinder.
- Clean tap water in a wash bottle.
- Liquid labeled T in a test tube.

#### **How to prepare liquid T.**

Dissolve reasonable amount of mixture of starch powder and glucose powder in warm water and then add liquid or live egg albumin and stir thoroughly to make a homogenous solution.

Teacher.co.ke