
Download this and other FREE materials from https://teacher.co.ke/notes

BIOLOGY (PRACTICALS) PAPER 3

231/3

CONFIDENTIAL



- Starch suspension labelled Liquid X
- Iodine solution
- Benedict's solution
- 2M hydrochloric acid (1ml)
- 2 Droppers
- Measuring cylinder (10 ml size)
- Means of heating/Bunsen burner
- 5 test-tubes
- Water in a small beaker
- Thermometer
- Test-tube holder
- 3 boiling tubes
- Tripod stand and gauze
- 3 labels
- White tile
- Water bath
- Diastase/ amylase enzyme (0.5g per student)

BN: -Liquid X is prepared by dissolving 5g of soluble starch in 50m1 of distilled water.

Thorough stirring is required whenever it is being used.

