

Name Adm No.

Class..... School.....

233/3

Date.....

CHEMISTRY

Paper 3

JUNE 2022

Time: 2¼ Hours

MARKING SCHEME

SUKELLEMO

FORM 4

CHEMISTRY PRACTICAL

INSTRUCTIONS

- * Write your **name** and the indicated details in the spaces provided on this page.
- * Answer **ALL** questions in the spaces provided in the question paper.
- * You are **NOT** allowed to start working with the apparatus for the *first 15 minutes* of the 2¼ hours allowed for this paper. This time is to enable you to read the question paper and make sure you have all the chemicals and apparatus you may need.
- * Mathematical tables and silent electronic calculators may be used.
- * All working must be clearly shown where necessary.

For Examiner's Use Only

Question	Maximum score	Candidate's score
1	21	
2	10	
3	9	
Total score	40	

This paper has 8 printed pages. Candidates should check the question paper to ascertain that ALL pages are printed as indicated and that no questions are missing.

QUESTION 1

You are provided with:

- 6.0g of an alkanolic acid labelled **solid P** in a boiling tube.
- 2M sodium hydroxide solution labelled **solution Q**.

You are required to:

- Determine the solubility of **solid P** at different temperatures.
- Determine the number of moles of water of crystallization in **solid P**.
- Find the molar mass of the alkanolic acid.

Procedure I

i) Using a burette add 10cm^3 of distilled water to the **solid P** in the boiling tube. Heat the mixture while stirring with a thermometer to about 70°C . When the entire solid has dissolved allow the solution to cool while stirring with the thermometer. Note the temperature at which crystals of **solid P** first appear. Record this temperature in the Table 1 below.

ii) Using the burette add 2cm^3 of distilled water to the contents of the boiling tube. Warm the mixture, while stirring with the thermometer until all the solid dissolves. Allow the mixture to cool while stirring. Note the temperature at which the crystals of **solid P** first appear.

iii) Repeat the procedure (ii) two more times and record the temperatures in Table I. **Retain the contents of the boiling tube for use in procedure II.**

TABLE I


Volume of water in the boiling tube (cm^3)	Temperature at which crystals of solid P first appear ($^\circ\text{C}$)	Solubility of Solid P ($\text{g}/100\text{gH}_2\text{O}$)
10	55	$\frac{6 \times 100}{10} = 60$ ✓
12	50	$\frac{6 \times 100}{12} = 50$ ✓
14	45	$\frac{6 \times 100}{14} = 42.9$ ✓
16	40	$\frac{6 \times 100}{16} = 37.5$ ✓

a) i) Complete Table 1 by calculating the solubility of **solid P** at different temperatures (6 mks)

- Complete table - 1mk
- Decimal point (whole number throughout / or 1dp) - 1mk
- Calculation - 2mk
- School value $\pm 2^\circ\text{C}$ - 1mk (whole school value is 1st reading) - 1mk
- Trend (decrease in temperature) - 1mk

6

ii) On the grid provided, plot a graph of the solubility of solid P (vertical axis) against temperature. (3 mks)



iii) Using the graph determine the temperature at which 55g of solid P would dissolve in 100cm³ of water. (1 mk)

52.5°C

Procedure II

- i) Transfer the contents of the boiling tube in Procedure I into a 250 ml volumetric flask. *Rinse both the boiling tube and the thermometer with distilled water and add to the volumetric flask.* Add more distilled water to make up to the mark. Transfer the solution to a 250 ml beaker, Label this solution R. Wash the volumetric flask then rinse it with distilled water ready for use in Step(ii).
- ii) Using a measuring cylinder place 25.0 cm³ of solution Q into a 250 ml volumetric flask. Add about 200cm³ of distilled water. Shake well. Add more distilled water to the make up to the mark. Label this solution T.
- iii) Fill a burette with solution R. Using a pipette and pipette filler, place 25 cm³ of solution T into a 250 ml conical flask. Add 2-3 drops of phenolphthalein indicator and titrate with solution R. Record your results in table II. Repeat the titration two more times and complete the table.

TABLE II

Titre	1	2	3
Final burette reading (cm ³)	13.2	27.2	13.2
Initial burette reading (cm ³)	0.0	14.0	0.0
Volume of solution R used (cm ³)	13.2	13.2	13.2

(4 marks)

CT-1
DP-1
ACC-1
PA-1
FA-1

DETERMINE:

Values averaged must be within ± 0.2 of each other

a) The average volume of solution R used (1 mk)

b) The concentration of solution T in moles per litre. (1 mk)

1000cm³ \Rightarrow 2mols

25cm³ $\Rightarrow \frac{25}{1000} \times 2 = 0.05 \text{ mols}$

0.05 mols \Rightarrow 250cm³

\Rightarrow 1000cm³

$\frac{1000}{250} \times 0.05 = 0.2 \text{ M}$

①

6

c) The concentration of alkanolic acid solution R in moles per litre (1 mole of acid reacts with 2 moles of the base) (2 mks)

<p>Moles of base in 25cm³ \rightarrow</p> $\frac{0.2 \times 25}{1000} = 0.005 \text{ moles}$	<p>Acid (2 mks)</p> $0.0025 \text{ moles} \Rightarrow \text{Av. volume}$	<p>Correct ans</p>
<p>Mole ratio of Acid: Base</p> $1 : 2$	<p>1000cm³</p> $= \frac{1000 \text{cm}^3}{\text{Av. volume}} \times 0.0025$	<p>Penalize -1/2 mk for wrong or missing units.</p>
<p>$\frac{0.005}{2} \times 1 = 0.0025 \text{ moles acid}$</p>		

d) The relative formula mass of the alkanolic acid, solid P (2 mks)

Ans in (c) above $\times 250$ \rightarrow correct ans.

$\frac{6}{\text{Ans in (c) above} \times 250} \times 1000$ OR $\frac{6}{\text{Ans above}}$ correct ans.

e) The formula of P has the form $M \cdot xH_2O$. Determine the value of x in the formula. (The relative formula mass of M is 90; O=16; H=1) (1mk)

RFM of P is $M \cdot xH_2O = \text{Ans in (d) above}$.

$x = \frac{\text{Ans in (d) above} - 90}{18}$

18

= correct ans.

Penalize -1/2 mk if the value of x is a fraction.

QUESTION 2

You are provided with **solid Y**. Use it to carry out the tests indicated below and record your observations and inferences in the spaces provided.

- a) Put half a spatula of **solid Y** in a boiling tube. Add about 5 cm³ of distilled water and shake. Divide the resultant mixture into four portions of 1 cm³ each.

Observations	Inferences
<p>1</p> <p>Solid dissolves to form a colourless solution (1 mk)</p>	<p>- Soluble salt (1 mk)</p> <p>Coloured ions absent (1 mk)</p> <p>Cu^{2+}, Fe^{2+}, Fe^{3+}</p>

- i) To the first portion add aqueous sodium hydroxide dropwise till in excess.

Observations	Inferences
<p>White precipitate dissolves in excess to form a colourless solution (1 mk)</p>	<p>Zn^{2+}, Al^{3+}, Pb^{2+} present (1 mk)</p>

- ii) to the second portion add aqueous ammonia dropwise till in excess

Observations	Inferences
<p>White precipitate dissolves in excess to form a colourless solution (1 mk)</p>	<p>Zn^{2+} present (1 mk)</p>

- iii) To the third portion add lead (II) nitrate solution and warm.

Observations	Inferences
<p>White precipitate does not dissolve on warming (1 mk)</p> <p>No yellow precipitate (1 mk)</p>	<p>SO_4^{2-}, SO_3^{2-}, CO_3^{2-} present (1 mk)</p> <p>Cl^- absent (1 mk)</p>

tied to correct observation

iv) To the fourth portion add dilute nitric (V) acid followed by Barium nitrate solution.

Observations	Inferences
- No effervescence ✓ (1/2) - white precipitate ✓ (1/2) (1 mk)	SO_4^{2-} present ✓ (1) (1 mk)

QUESTION 3

You are provided with solid M and solid F.

Carry out the tests below and write the observations and inference in the spaces below:

a) Using a metallic spatula, ignite about one half of solid M in a Bunsen burner flame.

Observations	Inferences
Solid melts and burns with a yellow sooty flame. (1) (1 mk)	$\text{C}=\text{C}$ present ✓ (1) $\text{C}\equiv\text{C}$ ✓ (1) (1 mk)

b) Place the other half of solid M in a boiling tube. Add 15cm³ of distilled water and shake well. Label this solution M. Use this solution for the following tests.

i) Place 2cm³ of solution M in test tube and determine pH.

Observations	Inferences
pH = 3 ✓ (1) Accept any value 1-3) not the range. (1/2 mk) (1/2 mk)	Strongly acidic ✓ (1) (1/2 mk)

ii) To about 2 cm³ of the solution M obtained in (b) above add 3 drops of acidified potassium manganate (VII)

Observations	Inferences
Purple KMnO_4 changes to colourless ✓ (1) (1 mk)	$\text{C}=\text{C}$; $\text{C}\equiv\text{C}$; ROH present ✓ (1/2) ✓ (1/2) (1 mk)

iii) To about 2 cm³ of the solution M obtained in (b) above add 2 drops of bromine water. *and warm*

Observations	Inferences
Yellow bromine water changes colourless ✓ (10) (1/2 mk)	>C=C< present ✓ (10) (1/2 mk)

iv) To the remaining solution M in the boiling tube, add half of solid ~~N~~ F

Observations	Inferences
Effervescence ✓ (10) (1/2 mk)	R-C(=O)-OH present ✓ (10) (1/2 mk) Accept RCOOH - COOH

2